




Expression microarray reagent guide

Reagents





Expression microarray reagent guide

Reagents for use with microarrays

GeneChip® 3' IVT PLUS Kit (IVT PLUS)^{1,2}

GeneChip® WT PLUS Reagent Kit (WT PLUS)¹

Affymetrix FlashTag™ Biotin HSR RNA Labeling Kit (FlashTag)

GeneChip® WT Pico Kit (WT Pico)

GeneChip® Pico Kit (Pico)¹

Arrays for use with WT PLUS¹

Clariom™ D arrays for human, mouse, and rat

Clariom™ S arrays for human, mouse, and rat

GeneChip® Human Transcriptome Array 2.0 (HTA 2.0)

GeneChip® human, mouse, and rat gene ST arrays

GeneChip® human, mouse, and rat exon 1.0 ST arrays

Arrays for use with WT Pico

Clariom™ D arrays for human, mouse, and rat

GeneChip® Human Transcriptome Array 2.0 (HTA 2.0)

GeneChip® human, mouse, and rat gene ST arrays²

GeneChip® human, mouse, and rat exon 1.0 ST arrays

Arrays for use with IVT PLUS^{1,3}

GeneChip® Human Genome U133 Plus 2.0 Array

GeneChip® PrimeView® Human Gene Expression Array

GeneChip® Mouse Genome 430 2.0 Array

GeneChip® Rat Genome 230 2.0 array

Arrays for use with Pico¹

Clariom™ S arrays for human, mouse, and rat

GeneChip® Human Genome U133 Plus 2.0 Array

GeneChip® PrimeView® Human Gene Expression Array

GeneChip® Mouse Genome 430 2.0 Array

GeneChip® Rat Genome 230 2.0 array

Arrays for use with FlashTag

GeneChip® miRNA arrays

List includes most common array types only and is not comprehensive. The most common arrays are paired with the appropriate reagent kits and sold as single assay part numbers for your convenience. Find additional assays, arrays, and information about the MyGeneChip™ Custom Array Program at www.affymetrix.com.

¹Also available as a high-throughput (HT) format for automation. ²Gene ST arrays should be paired with GeneChip® HT Pico Kit for automated sample preparation of challenging samples types and limited RNA input. ³Requires mRNA globin reduction for blood samples.

For Research Use Only. Not for use in diagnostic procedures.

P/N EMI06269 Rev. 2

© 2016 Affymetrix, Inc. All rights reserved. Affymetrix®, Axiom®, GeneChip®, Clariom™, CoMAP™, Command Console®, CytoScan®, DMET™, Eureka™, Eureka Genomics®, Eureka Genotyping™, Expression Console™, GeneAtlas®, GeneChip-compatible™, GeneTitan®, Genotyping Console™, myDesign™, MyGeneChip™, NetAffx®, OncoScan®, Powered by Affymetrix™, PrimeView®, and ViewRNA® are trademarks or registered trademarks of Affymetrix, Inc. Please see affymetrix.com/trademarks for a complete list of Affymetrix trademarks. All other trademarks are the property of their respective owners.